



Ready to use, aqueous solution for the rapid wiping disinfection of particularly sensitive surfaces and medical devices

Fields of application:

Acrylan® is a ready to use solution for the quick wiping disinfection of surfaces. The product distinguishes itself through its **excellent material compatibility**. Acrylan® is based on a QAC (quaternary ammonium compound) and is free of alcohols and aldehydes. Therefore, it is ideally suited for the application onto alcohol-sensitive materials like acrylic glass.

Acrylan® offers a broad spectrum of efficacy and good cleaning properties. It is effective within short contact times and has a pleasant, neutral scent.

The product is suitable for the disinfection and cleaning of medical devices and surfaces, e. g. in ICUs, examination tables, stretchers, and protective mattress covers in doctors' practices, accident and emergency units and physiotherapy; also for ambulance surfaces, ultrasound devices, tonometer measuring prisms and contact glasses (e. g., by Haag-Streit), etc.

Recommendations for use:

Test	Concentration	Contact time
VAH – Surface Disinfection with mechanical action I VAH Standard Methods		
VAH – Surface Disinfection: bactericidal, yeasticidal clean and dirty conditions (equivalent to EN 13727, EN 13624 and EN 16615)	undiluted	5 min
DVV / RKI Guideline – virucidal efficacy		
limited virucidal* (inactivation of enveloped viruses) (incl. HIV / HBV / HCV)	undiluted	5 min
Vacciniavirus	undiluted	5 min
BVDV (Bovine Viral Diarrhea Virus)	undiluted	1 min
Rotavirus	undiluted	5 min

* according to DVV/RKI recommendation, Bundesgesundheitsblatt 03/2017

Special characteristics:

- suitable for the use in the wipes dispensing system **A°-Wipes** and **Kimtech® Wettask® DS**
- wide range of application
- fast onset of action
- very high material compatibility even with particularly sensitive materials like acrylic glass
- pleasant scent

Composition (active ingredients):

100 g contain:
0.525 g Didecyltrimethylammonium chloride
0.120 g Polyhexanide

Listing:

- VAH-list
- ÖGHMP-list

Spectrum of efficacy:

- bactericidal
- yeasticidal
- limited virucidal*

Antisepsis by Antiseptica – Antisepsis mit Vernunft

The company is certified according to **EN ISO 13485** and meets the Council's requirements of Regulation (EU) 2017/745 of the European Parliament and the Council on Medical Devices.

Lysoform Dr. Hans Rosemann GmbH, Kaiser-Wilhelm-Straße 133, D-12247 Berlin, Telefon: +49 30 77992-0, Fax: +49 30 77992-219, www.lysoform.de



ANTISEPTICA Dr. Hans-Joachim Molitor GmbH
D-12247 Berlin, Kaiser-Wilhelm-Str. 133
Tel.: +49 30 77992-200
Mail: info@antiseptica.com, www.antiseptica.com

Distribution in Austria by:
ANTISEPTICA chem.-pharm. Produkte GmbH
A-1090 Wien, Frankgasse 6 / Top 5
Tel.: +43 1 3746600, Fax: +43 1 374 660066
Mail: office@antiseptica.at, www.antiseptica.at

Physical and chemical properties:

Appearance: clear, colourless liquid
 Odour: neutral
 pH-value: approx. 7 at 20 °C
 Density: approx. 1 g/cm³
 Flash point: > 60 °C according to DIN 51755

Application:

Surface disinfection with mechanical action;
 scrub-wipe disinfection:

Apply a sufficient amount (wet-wiping) of the undiluted, ready to use product to the surface by using a scrub-wipe disinfection with mechanical distribution. Ensure a thorough, complete moistening of the surface and allow to operate for the recommended contact time.

Massive contamination with organic material, (e. g., blood, secretions), should be removed before disinfection using a disposable cloth moistened with disinfectant, which should then be properly disposed.

For use in combination with the wipes dispensing systems A®-Wipes: Acrylan® can be used for a shelflife of 28 days.

Special advice:

- Wear adequate personal protective equipment when handling disinfectants.
- Avoid contact with eyes and mucous membranes.
- Clean and rinse surfaces, buckets and wiping utensils properly before first use or when switching products.
- Do not mix with other substances (e. g. anionic surfactants, aldehydes).
- Please follow local guidelines for the cleaning and disinfection of surfaces!
- Ensure adequate ventilation.
- Protect from sunlight.

Forms of delivery:

- carton of 10 x 1 L plastic bottles
- carton of 2 x 5 L plastic canisters

also available:

- various application techniques that can be screwed onto bottles and canisters

When handling the product, please note the information according to the Hazardous Substances Ordinance given in the safety data sheet.

Further information is provided by your expert advisor upon request or on our website.

Use disinfectants safely. Always read the label and product information before use.

Note for disposal:

The solution of Acrylan® can be disposed of with waste water, as it is common practice for cleaning agents. A neutralisation or inactivation before release into sewage is neither necessary nor ecologically worthwhile.

CE 0482

Antisepsis by Antiseptica – Antisepsis mit Vernunft

The company is certified according to **EN ISO 13485** and meets the Council's requirements of Regulation (EU) 2017/745 of the European Parliament and the Council on Medical Devices.

Lysoform Dr. Hans Rosemann GmbH, Kaiser-Wilhelm-Straße 133, D-12247 Berlin, Telefon: +49 30 77992-0, Fax: +49 30 77992-219, www.lysoform.de



ANTISEPTICA Dr. Hans-Joachim Molitor GmbH
 D-12247 Berlin, Kaiser-Wilhelm-Str. 133
 Tel.: +49 30 77992-200
 Mail: info@antiseptica.com, www.antiseptica.com

Distribution in Austria by:
ANTISEPTICA chem.-pharm. Produkte GmbH
 A-1090 Wien, Frankgasse 6 / Top 5
 Tel.: +43 1 3746600, Fax: +43 1 374 660066
 Mail: office@antiseptica.at, www.antiseptica.at