

Antiseptica Combi Spray

ANTISEPTICA
Hygiene-Lösungen seit 1976
Solutions for Hygiene since 1976



Alcoholic solution
for the targeted, rapid disinfection of surfaces
with very broad efficacy in very short contact times

FIELDS OF APPLICATION:

It is well known that direct contact to contaminated surfaces may spread pathogens. Therefore, a consistent surface disinfection is indispensable, also under great time pressure. In this case, alcoholic rapid disinfectants, which are effective in no more than one minute, are often selected.

In addition, the demand for active ingredients that affect staff, patients and the environment as little as possible while being reliably effective remains.

Antiseptica Combi Spray is a ready to use, alcoholic solution for the reliable **rapid surface disinfection** with a broad spectrum of efficacy even against non-enveloped viruses.

The formaldehyde-free product is suitable for the disinfection of medical inventory and medical devices like operating inventory and operating tables, dressing trolleys, transport carts and trays, as well as in important functional areas such as operating theatres, ICUs, and ambulances, in doctors' and dentists' practices as well as in secondary areas of human medicine.

Antiseptica Combi Spray is also suitable for disinfection measures in the pharmaceutical and cosmetic industry, as well as in laboratories and diagnostics.

SPECIAL CHARACTERISTICS:

- very broad spectrum of efficacy
- rapid onset of action
- dries quickly
- pleasant scent
- free of formaldehyde

RECOMMENDATIONS FOR USE:

Test	Concentration	Contact time
VAH – Surface disinfection with mechanical action VAH Standard Methods		
VAH – Surface disinfection: bactericidal, yeasticidal clean and dirty conditions (equivalent to EN 13727 and EN 13624)	undiluted	5 min
Surface disinfection: bactericidal, yeasticidal clean and dirty conditions	undiluted	1 min
Surface disinfection with mechanical action DGHM Standard Methods		
tuberculocidal clean and dirty conditions (equivalent to EN 14348)	undiluted	1 min
EN 14476 – virucidal efficacy		
virucidal activity against enveloped viruses	undiluted	5 min
limited spectrum of virucidal activity	undiluted	5 min
virucidal	undiluted	60 min

*according to RKI-recommendation, Bundesgesundheitsblatt 3/2017

Antisepsis by Antiseptica – Antisepsis mit Vernunft

The company is certified according to **EN ISO 13485** and meets the Council's requirements of Regulation (EU) 2017/745 of the European Parliament and the Council on Medical Devices.

 **Lysoform Dr. Hans Rosemann GmbH**, Kaiser-Wilhelm-Straße 133, D-12247 Berlin, Telefon: +49 30 77992-0, Fax: +49 30 77992-219, www.lysoform.de



 **ANTISEPTICA Dr. Hans-Joachim Molitor GmbH**
D-12247 Berlin, Kaiser-Wilhelm-Str. 133
Tel.: +49 30 77992-200
Mail: info@antiseptica.com, www.antiseptica.com

Distribution in Austria by:
ANTISEPTICA chem.-pharm. Produkte GmbH
A-1090 Wien, Frankgasse 6 / Top 5
Tel.: +43 1 3746600, Fax: +43 1 374 660066
Mail: office@antiseptica.at, www.antiseptica.at

APPLICATION:

Surface disinfection with mechanical action; scrub-wipe disinfection:

Apply a sufficient amount (wet-wiping) of the undiluted, ready to use product to the surface by using a scrub-wipe disinfection with mechanical distribution. Ensure a thorough, complete moistening of the surface and allow to operate for the recommended contact time.

Massive contamination with organic material (e. g. blood, secretions) should be removed before disinfection using a disposable cloth moistened with disinfectant, which should then be properly disposed.

PHYSICAL AND CHEMICAL PROPERTIES:

Appearance:	clear, colourless liquid
Odour:	alcoholic
pH-value:	approx. 6 at 20 °C
Density:	approx. 0.9 g/cm ³ at 20 °C
Flash point:	28 °C according to DIN 51755

SPECTRUM OF EFFICACY:

- bactericidal
- yeasticidal
- tuberculocidal
- **virucidal***

LISTING:

- VAH-list
- ÖGHMP-list

COMPOSITION (ACTIVE INGREDIENTS):

100,00 g contain:
35.00 g Propan-1-ol
20.00 g Ethanol
0.05 g Glutaraldehyde

SPECIAL ADVICE:

- Alcoholic disinfectants should not be used for the disinfection of large surfaces (danger of fire and explosion!). Output: max. 50 ml/m².
- Keep away from sources of ignition; no smoking. Hot surfaces must be cooled down before disinfection. Electric equipment must be switched to neutral/off and volt-free.
- Contains volatile substances. Ensure adequate ventilation.
- Not suitable for the disinfection of alcohol-sensitive surfaces (e. g. acrylic glass, alcohol-soluble lacquers).
- Avoid contact with eyes, skin, and mucous membranes. Do not breathe vapour.
- Wear adequate personal protective equipment when handling disinfectants.
- Clean and rinse surfaces, buckets and wiping utensils properly before first use or when switching products.
- Do not mix with other substances (e. g. strong bases, oxidising agents).
- Please follow local guidelines for the cleaning and disinfection of surfaces!

FORMS OF DELIVERY:

- carton of 10 x 1 L plastic bottles
- carton of 2 x 5 L plastic canister
- 1 x 10 L canister

also available:

- various application techniques that can be screwed onto bottles and canisters

Further information is provided by your expert advisor upon request or on our website.

Use disinfectants safely. Always read the label and product information before use.


Surface disinfectant for medical devices according to Regulation (EU) 2017/745 „Medical Device Regulation (MDR)“.

CE 0482

Biocidal product according to BPR
BAuA Reg.-No.: N-45588

Antisepsis by Antiseptica – Antisepsis mit Vernunft

The company is certified according to **EN ISO 13485** and meets the Council's requirements of Regulation (EU) 2017/745 of the European Parliament and the Council on Medical Devices.

 **Lysoform Dr. Hans Rosemann GmbH**, Kaiser-Wilhelm-Straße 133, D-12247 Berlin, Telefon: +49 30 77992-0, Fax: +49 30 77992-219, www.lysoform.de



 **ANTISEPTICA Dr. Hans-Joachim Molitor GmbH**
D-12247 Berlin, Kaiser-Wilhelm-Str. 133
Tel.: +49 30 77992-200
Mail: info@antiseptica.com, www.antiseptica.com

Distribution in Austria by:
ANTISEPTICA chem.-pharm. Produkte GmbH
A-1090 Wien, Frankgasse 6 / Top 5
Tel.: +43 1 3746600, Fax: +43 1 374 660066
Mail: office@antiseptica.at, www.antiseptica.at