

Dosage table for liquids

	0.25 %	0.50 %	1.00 %	1.50 %	2.00 %	2.50 %	3.00 %	4.00 %
1 litre	2.5 ml	5 ml	10 ml	15 ml	20 ml	25 ml	30 ml	40 ml
2 litres	5 ml	10 ml	20 ml	30 ml	40 ml	50 ml	60 ml	80 ml
3 litres	7.5 ml	15 ml	30 ml	45 ml	60 ml	75 ml	90 ml	120 ml
4 litres	10 ml	20 ml	40 ml	60 ml	80 ml	100 ml	120 ml	160 ml
5 litres	12.5 ml	25 ml	50 ml	75 ml	100 ml	125 ml	150 ml	200 ml
6 litres	15 ml	30 ml	60 ml	90 ml	120 ml	150 ml	180 ml	240 ml
7 litres	17.5 ml	35 ml	70 ml	105 ml	140 ml	175 ml	210 ml	280 ml
8 litres	20 ml	40 ml	80 ml	120 ml	160 ml	200 ml	240 ml	320 ml
9 litres	22.5 ml	45 ml	90 ml	135 ml	180 ml	225 ml	270 ml	360 ml
10 litres	25 ml	50 ml	100 ml	150 ml	200 ml	250 ml	300 ml	400 ml
15 litres	37.5 ml	75 ml	150 ml	225 ml	300 ml	375 ml	450 ml	600 ml
20 litres	50 ml	100 ml	200 ml	300 ml	400 ml	500 ml	600 ml	800 ml
25 litres	62.5 ml	125 ml	250 ml	375 ml	500 ml	625 ml	750 ml	1,000 ml
30 litres	75 ml	150 ml	300 ml	450 ml	600 ml	750 ml	900 ml	1,200 ml
40 litres	100 ml	200 ml	400 ml	600 ml	800 ml	1,000 ml	1,200 ml	1,600 ml
50 litres	125 ml	250 ml	500 ml	750 ml	1,000 ml	1,250 ml	1,500 ml	2,000 ml

Instructions:

First pour the water at approx. room temperature into the bucket and then add the disinfectant. In order to obtain the correct volume of water to be added, subtract the volume of the required disinfectant from the volume of the desired dilute solution (first column), e.g. 50 litres of a 4 % dilute solution => 48 litres water + 2 litres (2,000 ml) disinfectant.



ANTISEPTICA

Carl-Friedrich-Gauß-Str.7, D-50259 Pulheim
 Tel. +43 - 2234 - 9 84 66 - 0, Fax +43 - 2234 - 9 84 66 - 11
 Internet: <http://www.antiseptica.com> E-Mail: info@antiseptica.com