



Under application for union authorisation  
according to regulation (EU) 528/2012 (BPR)!

Ready to use, unperfumed alcoholic  
hand rub for hygienic and surgical  
hand disinfection

#### Fields of application:

Nosocomial infections are most commonly transmitted by contaminated hands. Thus, regular and thorough hand disinfection is vital for effective infection prevention and control, particularly in the medical sector. In order to improve compliance in hand hygiene, a hand disinfectant should efficiently reduce pathogens in a short contact time, be effective quickly and provide very good skin compatibility.

Manorapid® pure is a **ready to use, alcoholic and perfume free hand rub**, for hygienic and surgical hand disinfection. The product possesses a strong microbicidal efficacy, whereas the caring agent lanolin provides **excellent skin compatibility**, making the product highly suitable for frequent use. The absence of perfume is beneficial to people with sensitive skin.

#### Special characteristics:

- Free from fragrances and colourants
- Moisturising
- Broad spectrum of activity
- Dries quickly
- Developed for professional users

#### Spectrum of efficacy:

- bactericidal
- yeasticidal
- tuberculocidal
- limited virucidal\*
- effective against Noro- and Rotavirus

#### Recommendations for use:

Test	Concentration	Contact time
<b>EN Standard Methods hand disinfection</b>		
hygienic hand disinfection (EN 1500)	undiluted	30 s
surgical hand disinfection (EN 12791)	undiluted	1,5 min.
<b>DVV / RKI Guideline – virucidal efficacy</b>		
limited virucidal* (inactivation of enveloped viruses) (incl. HIV / HBV / HCV)	undiluted	30 s
Vaccinia virus	undiluted	30 s
BVDV (Bovine Viral Diarrhea Virus)	undiluted	30 s
Coronavirus (SARS)	undiluted	30 s
Adenovirus	undiluted	2 min.
<b>EN 14476 – virucidal efficacy</b>		
limited virucidal* (clean conditions)	undiluted	30 s
Norovirus (MNV) <sup>a</sup> (clean conditions)	undiluted	30 s
Rotavirus (clean conditions)	undiluted	30 s
<b>EN 14348 – tuberculocidal efficacy</b>		
tuberculocidal (clean conditions)	undiluted	30 s

\* according to DVV/RKI declaration, Bundesgesundheitsblatt 3/2017, <sup>a</sup> murine Norovirus = surrogate virus for human Norovirus

## Antiseptics by Antiseptica - Made in Germany

The company is certified according to **EN ISO 13485** and meets the Council's requirements of Regulation **(EU) 2017/745** of the European Parliament and the Council on Medical Devices.



■ Lysoform Dr. Hans Rosemann GmbH

Kaiser-Wilhelm-Straße 133, D-12247 Berlin, phone +49 30 77992-0, fax +49 30 77992-219, [www.lysoform.de](http://www.lysoform.de)

ANTISEPTICA Dr. Hans-Joachim Molitor GmbH

Kaiser-Wilhelm-Str. 133, D-12247 Berlin, phone +49 30 77992-200, mail: [info@antiseptica.com](mailto:info@antiseptica.com), [www.antiseptica.com](http://www.antiseptica.com)

**Application:****For hygienic hand disinfection:**

Rub undiluted product into hands and keep them moist for the entire contact time.

**For surgical hand disinfection:**

Rub undiluted product into hands and lower arms and keep them moist for the entire contact time.

**Physical and chemical characteristics:**

Appearance: clear, colourless liquid

Scent: alcoholic

pH-value: 6 – 8 at 20 °C

Density: approx. 0.9 g/cm<sup>3</sup> at 20 °C

Flash point: 19.5 °C according to DIN 51755

**Composition (active ingredients):**

100 ml contain: 63.1 g Propan-2-ol

**Forms of delivery:**

- carton of 20 x 500 ml bottles
- carton of 12 x 1 L bottles
- carton of 2 x 5 L canister

also available:

- wall dispenser for 500 ml bottles
- wall dispenser for 1 L bottles
- dosing pump for 500 ml bottles
- dosing pump for 1 L bottles
- dosing pump for 5 L canister

When handling the product, please note the information according to the Hazardous Substances Ordinance given in the safety data sheet.

Further information is provided by your expert advisor upon request or on our website.

Use disinfectants safely. Always read the label and product information before use.

Biocidal product according to BPR  
BAuA Reg.-No.: N-65660

## Antiseptica by Antiseptica - Made in Germany

The company is certified according to **EN ISO 13485** and meets the Council's requirements of Regulation **(EU) 2017/745** of the European Parliament and the Council on Medical Devices.



■ Lysoform Dr. Hans Rosemann GmbH

Kaiser-Wilhelm-Straße 133, D-12247 Berlin, phone +49 30 77992-0, fax +49 30 77992-219, [www.lysoform.de](http://www.lysoform.de)

ANTISEPTICA Dr. Hans-Joachim Molitor GmbH

Kaiser-Wilhelm-Str. 133, D-12247 Berlin, phone +49 30 77992-200, mail: [info@antiseptica.com](mailto:info@antiseptica.com), [www.antiseptica.com](http://www.antiseptica.com)